



## Pure DNA, Pure Data


### AGENCOURT® AMPURE®

PCR PURIFICATION SYSTEM

Agencourt AMPure is a highly efficient PCR\* purification system that delivers superior quality DNA. The Agencourt AMPure method utilizes Solid Phase Reversible Immobilization (SPRI®) magnetic bead-based technology which requires no centrifugation or filtration. It can be easily used in manual and automated 96- and 384-well formats. Using a simple, three-step protocol, reaction contaminants are removed to improve downstream performance in applications such as sequencing and SNP genotyping.

Table 1

	Agencourt AMPure	ExoSAP-IT
Number of Samples	16	16
Pass Rate**	100%	94%
Phred20	701±45	542±142



**Figure 1.** A 967 bp fragment of a pGEM plasmid was amplified, purified using Agencourt AMPure (A) or ExoSAP-IT (B) and sequenced using BigDye Terminator v3.1. Unincorporated dye terminators were removed using Agencourt CleanSEQ® and samples detected on an ABI PRISM 3730xl.

### Efficient PCR Purification

Agencourt AMPure provides:

- Efficient removal of unincorporated dNTPs, primers, salts and other contaminants
- High recovery of both small and large amplicons
- A scalable and automation compatible method

### Flexible Process

Agencourt AMPure flexible protocol is compatible with both 96- and 384-well formats and can be used for a variety of applications:

- PCR Sequencing
- SNP Genotyping
- cDNA Labeling
- Rolling Circle Amplification


### High Purity – High Performance

The Agencourt AMPure purification process binds PCR amplicons to para-magnetic particles and draws them out of solution, allowing contaminants such as primers, salts, and dNTPs to be easily rinsed away. The purified PCR product is an ideal template for sequencing. Alternative PCR treatments, such as ExoSAP-IT, merely denature enzymes and degrade unincorporated primers and dNTPs without removing them. 976bp PCR amplicons were generated from pGEM plasmid and Qiagen HotStar PCR mix, and the resulting PCR reactions were purified using either Agencourt AMPure or treated with ExoSAP-IT. Table 1 shows that Agencourt AMPure gives longer phred20 readlengths and lower standard deviation than ExoSAP-IT. The poor purity of the ExoSAP-IT template interferes with the sequencing reaction resulting in shorter readlengths, decreased quality and increased background. Figure 1 shows Agencourt AMPure traces versus ExoSAP-IT traces. Agencourt AMPure's efficient removal of contaminants results in higher quality sequencing data.


### Effective Cleanup and Reliable Results

The Agencourt AMPure purification system can be used to purify amplicons from 100bp to 40kb. This method routinely produces recovery rates greater than 90% for





**Figure 2.** 165, 321, 410, 508, and 967 bp fragments of the amplified pGEM plasmid. Equal amounts of samples were run on a 4% agarose gel. Lane MW: 100 bp ladder; Lanes 1, 3, 5, 7, and 9: unpurified PCR products; Lanes 2, 4, 6, 8, and 10: Agencourt AMPure-purified samples.



**Figure 3.** Eight 967 bp (A) and 165 bp (B) PCR products were amplified. Samples were purified using Agencourt AMPure and separated by electrophoresis on a 1.2% agarose gel. Lanes MW: 100 bp MW ladder; Lane 1: Unpurified PCR products; Lanes 2 - 9: Agencourt AMPure purified samples.

amplicons larger than 200 bp and greater than 85% for smaller PCR amplicons (Figure 2). The utilization of a magnetic bead-based method allows Agencourt® AMPure® to generate reproducible results when purifying large and small amplicons. Figure 3 shows two gels following the Agencourt AMPure purification of 967 bp and 165 bp PCR fragments demonstrating consistent recovery of both small and large amplicons.

### Agencourt AMPure Automation

The use of magnetic beads eliminates centrifugation and filtration allowing Agencourt AMPure protocols to be completely automated while maintaining superior product quality. Agencourt AMPure is fully compatible with several common automation liquid handling platforms, and is optimized for Beckman Coulter's Biomek® series of automated laboratory workstations. Contact Agencourt for the most up-to-date list of supported platforms.



### Ordering Information

For product pricing, please visit our website at [www.agencourt.com](http://www.agencourt.com) or contact your local sales representative.

Product	Size***	Product #
Agencourt AMPure 96 Starter Kit	3,333 preps	000146
Agencourt AMPure 384 Starter Kit	3,333 preps	000148
Agencourt AMPure 60 mL Kit	3,333 preps	000130
Agencourt AMPure 450 mL Kit	25,000 preps	000132
Related Products		Product #
Agencourt SPRIPlate® 96R Magnet Plate		000219
Agencourt SPRIPlate 384 Magnet Plate		000222

\* The PCR process is covered by patents owned by Roche Molecular Systems, Inc., and F Hoffman-La Roche, Ltd.  
 \*\* Passing Read: Average P<sub>read</sub> value greater than 20 for bases between 100bp and 300bp.  
 \*\*\* Typical reaction volume is 10–20 µL for the 96-well microplate and 5–10 µL for the 384-well microplate. Number of preps based on 10 µL reaction volume.  
 All trademarks are the property of their respective owners.



Innovate Automate  
SIMPLIFY

Agencourt Bioscience Corporation, A Beckman Coulter Company • 800-361-7780 • [www.agencourt.com](http://www.agencourt.com)  
 500 Cummings Center, Suite 2450 • Beverly, Massachusetts 01915

#### Worldwide Offices:

Australia (61) 2 9844-6000 Canada (905) 819-1234 China (86) 10 6515 6028 Eastern Europe, Middle East, North Africa (41) 22 994 07 07  
 France 01 49 90 90 00 Germany (31) 10 470 79 26 Hong Kong (852) 2814 7431/2814 0481 Italy 02-953921 Japan 03-5404-8359  
 Mexico (55) 560-57770 Netherlands (31) 10 470 79 26 Singapore (65) 6339 3633 South Africa, Sub-Saharan Africa (27) 11-805-2014/5 Spain 91 3836080  
 Sweden 08-564 85 900 Switzerland 0800 850 810 Taiwan (886) 2 2378 3456 Turkey 90 216 309 1900 U.K. 01494 441181 U.S.A. 1-978-867-2600