



High-Performing DNA Centrifugation System

ALLEGRA™ 25R BENCH CENTRIFUGE

WITH NEW S5700 AND TA-15-1.5 DNA ROTORS

For a quick and easy solution for DNA isolation/purification and post-reaction cleanup, use the Allegra 25R General Purpose Centrifuge with S5700 and TA-15-1.5 rotors.

DNA ISOLATION AND PURIFICATION

- Manual DNA mini-preps or spin filtration columns – up to thirty 1.5 mL tubes at $25,160 \times g$
- Standard 96-well plates, deep-well plates and square-well plates, DNA kit 96-well filtration plates – up to $6,130 \times g$

POST-REACTION CLEANUP

- PCR plates, DNA kit 96-well filtration plates, deep-well and square-well plates – up to $6,130 \times g$
- Plate Precipitations – **7.5 MINUTES** total centrifuge time – start to finish

TOTAL COMPLEMENT OF ROTORS


- Pellet bacteria, yeast cell cultures or proteins
 - 4×500 mL – up to $5,530 \times g$
 - 6×250 mL – up to $15,300 \times g$
 - 8×50 mL – up to $21,000 \times g$



*TA-15-1.5
DNA Rotor*


S5700 DNA Rotor

FOR USE IN THE NEW S5700 DNA ROTOR




Standard 96-Well Plates

Per adapter = 5
Per rotor = 10




Deep-Well Plates

Per adapter = 2
Per rotor = 4




Square-Well Plates

Per adapter = 1
Per rotor = 2



DNA Kit 96-Well Filtration Plates^{††}

Per adapter = 1
Per rotor = 2



PCR Plates^{†††}

Per adapter = 1
Per rotor = 2

ORDERING INFORMATION

Allegra™ 25R Centrifuges

PART NO.	DESCRIPTION
369434	Allegra 25R High-Speed Benchtop Centrifuge (60 Hz, 208V)
369435	Allegra 25R High-Speed Benchtop Centrifuge (50/60 Hz, 200V)
369436	Allegra 25R High-Speed Benchtop Centrifuge (50 Hz, 230V)

Rotors

PART NO.	DESCRIPTION	CAPACITY	RPM	g-FORCE
368954	S5700 Swinging-Bucket Rotor (Includes part nos: 368920, 368957, 267002, 267005)	2 × 5 Std., 2 × 2 Deep-Well, 2 × 1 Square-Well, 2 × 1 PCR, 2 × 1 96-Well Filtration Plates	5,700	6,130
368298	TA-15-1.5 Fixed-Angle Rotor	30 × 0.25-2.0 mL	15,000	25,160
368308	TS-5.1-500 Swinging-Bucket Rotor	4 × 500 mL	5,100	5,530
368293	TA-10-250 Fixed-Angle Rotor	6 × 250 mL	10,000	15,300
368303	TA-14-50 Fixed-Angle Rotor	8 × 50 mL	14,500	21,100

S5700 Adapters

PART NO.	DESCRIPTION
368920	Microplate Carrier Kit (Set of 2)
368957	Support Pad (Set of 4)
267002	96-Well Cap Mat (Nonsterile)
267005	96-Well Cap Mat (Sterile)

TS-5.1-500 Adapters

PART NO.	DESCRIPTION	CONFIGURATION
368451	Microplate Carriers (Set of 2), 4,100 rpm; 3,010 × g	4 × 2 Std., 4 × 1 Deep-Well Plate
368464	15 mL Falcon Tube (Conical/Set of 2)	9 × 15 mL, 70 × 120 mm
368461	50 mL Falcon Tube (Conical/Set of 2)	4 × 50 mL, 30 × 120 mm

TA-15-1.5 Adapters

PART NO.	DESCRIPTION	CONFIGURATION
364701	1.5 mL Polyallomer Tube (Set of 12)	1 × 1.5 mL, 11 × 40 mm
361247	0.4 mL Polyethylene Tube (24/pk)	1 × 0.4 mL, 7 × 40 mm
361247	0.25 mL Polyethylene Tube	1 × 0.25 mL, 7 × 30 mm

For detailed rotor information, visit us at
www.beckmancoulter.com/centrifugefirst



I.S. EN ISO 9001

Developing innovative solutions in genetic analysis, drug discovery, and instrument systems



Innovate SIMPLIFY **Automate**

Beckman Coulter, Inc. • 4300 N. Harbor Boulevard, Box 3100 • Fullerton, California 92834-3100
Sales: 1-800-742-2345 • Service: 1-800-551-1150 • Telex: 678413 • Fax: 1-800-643-4366 • www.beckmancoulter.com

Worldwide Bioresearch Division Offices:

Australia (61) 2 9844-6000 **Canada** (905) 819-1234 **China** (86) 10 6515 6028 **Eastern Europe, Middle East, Africa** (41) 22 994 07 07
France 01 49 90 90 00 **Germany** (89) 35870-0 **Hong Kong** (852) 2814 7431/2814 0481 **Italy** 02-953921 **Japan** 03-5404-8359
Mexico 525-559-16-35 **Netherlands** 0297-230630 **Singapore** (65) 339 3633 **South Africa** (27) 11-805-2014/5 **Spain** 91 728 7900
Sweden 08-564 85 900 **Switzerland** 0800 850 810 **Taiwan** (886) 2 2378 3456 **U.K.** 01494 441181 **U.S.A.** 1-800-742-2345.